

EXTREMITY WHIRLPOOLS

A E SERIES

- 1** • Dimensions: 25" x 13" x 15"

15 Gallon **Description**

070513 Mobile

070512 Mobile With Sump

070514 Stationary

- 2** • Dimensions: 28" x 15" x 18"

22 Gallon **Description**

070516 Mobile

070515 Mobile With Sump

070517 Stationary

- 3** • Dimensions: 28" x 15" x 21"

27 Gallon **Description**

070519 Mobile

070518 Mobile With Sump

070520 Stationary

- 4** • Dimensions: 28" x 15" x 26"

36 Gallon **Description**

070522 Mobile

070521 Mobile With Sump

070523 Stationary

- 5** • Dimensions: 32" x 15" x 25"

45 Gallon **Description**

070525 Mobile

070524 Mobile With Sump

070526 Stationary

B FERNO-ILLE WHIRLPOOLS

- 6** • Suggested use for foot or ankle.
• Dimensions: 24" L x 14" W x 14" D

Description

070051 14 Gallon; Mobile

070052 14 Gallon; Stationary

- 7** • Suggested use for foot or ankle.
• Dimensions: 24" L x 14" W x 16" D

Description

070053 16 Gallon; Mobile

070054 16 Gallon; Stationary

- 8** • Suggested use for lower leg, knee or arm.
• Dimensions: 28" L x 16" W x 25" D

Description

070057 36 Gallon; Mobile

070058 36 Gallon; Stationary

- 9** • These tanks are manufactured at an intermediate height in order to accommodate both upper and lower extremities.

- Suggested use for foot, ankle or arm.
• Dimensions: 28"L x 16"W x 18"D

Description

070055 22 Gallon; Mobile

070056 25 Gallon; Stationary

C LI'L CHAMP AQUA WHIRL

Description

070500 4.5 Gallon; 23.25" x 13.75" x 19.25"

070502 Replacement Turbine

070501 4.5 Gallon; with Separate Drain Pump



Ferno-Ille mobile tanks can be easily rolled to water source or treatment area and then later stored out of the way. No permanent plumbing connections required.

Tanks can be completely drained with a powered drain pump into any accessible drain.



Li'l Champ is a compact whirlpool tank designed for podiatry and hand-therapy use on table-tops or floors.



ADJUSTABLE HEIGHT WHIRLPOOLS

A FERNO-ILLE WHIRLPOOLS

- 1** • Dimensions: 24”L x 14”W x 14”D

Description

070059 14 Gallon; Stainless-Steel

- 2** • Approximate height range: 14” adjustable.
• Suggested use for foot, ankle or arm.
• Dimensions: 38”L x 23”W x 13.5”D

Description

070060 12 Gallon; Fiberglass

B AQUA WHIRL WHIRLPOOLS

- 3** • Dimensions: 30” x 20” x 14.5”

Description

070570 Plastic Tank

- 4** • Dimensions: 24” x 14.5” x 15”

Description

070571 Stainless Steel Tank

HIGH FULL BODY WHIRLPOOLS

C FERNO-ILLE WHIRLPOOLS

- Offers complete coverage to different areas of the body.
- If full immersion is not desired, an optional adjustable suspension seat can be utilized.
- Suggested use for lower back, leg and hip.
- Dimensions: 36”L x 22”W x 28”D

Description

070061 72 Gallon; Mobile

070063 110 Gallon; Mobile

070062 72 Gallon; Stationary

070064 110 Gallon; Stationary

D H SERIES HI-BOY

- 5** • Dimensions: 36” x 20” x 28”

60 Gallon

Description

070535 Mobile

070534 Mobile with Sump

070536 Stationary

- 6** • Dimensions: 42” x 20” x 28”

75 Gallon

Description

070538 Mobile

070537 Mobile with Sump

070539 Stationary

- 7** • Dimensions: 48” x 20” x 28”

90 Gallon

Description

070541 Mobile

070540 Mobile with Sump

070542 Stationary

- 8** • Dimensions: 48” x 24” x 28”

105 Gallon

Description

070544 Mobile

070543 Mobile with Sump

070545 Stationary



070553

E



F



These long and shallow units are designed for easy access.

070562

G



070563



LONG AND LOW FULL BODY WHIRLPOOLS

E L SERIES LO-BOY

- 9** • Dimensions: 52" x 24" x 18"

75 Gallon

Description

070547 Mobile

070546 Mobile with Sump

070548 Stationary

- 10** • Dimensions: 60" x 24" x 18"

90 Gallon

Description

070550 Mobile

070549 Mobile with Sump

070551 Stationary

- 11** • Dimensions: 66" x 24" x 18"

105 Gallon

Description

070553 Mobile

070552 Mobile with Sump

070554 Stationary

F FERNO-ILLE WHIRLPOOLS

- 12** • The lower height makes these ideal for pediatric and geriatric patients requiring hydrotherapy treatment.
- Suggested use for upper back, leg and hip.
 - Dimensions: 66"L x 25"W x 18"D

Description

070067 102 Gallon; Mobile

070068 102 Gallon; Stationary

- 13** • The lower height makes these ideal for pediatric and geriatric patients requiring hydrotherapy treatment.
- Suggested use for upper back, leg and hip.
 - Dimensions: 60"L x 25"W x 18"D

Description

070065 93 Gallon; Mobile

070066 93 Gallon; Stationary

SPORTS TRAINER WHIRLPOOLS

G S SERIES

- 14** • Dimensions: 56" x 24" x 25"

110 Gallon

Description

070562 Mobile

070561 Mobile with Sump

070563 Stationary

- 15** • Dimensions: 48" x 20" x 25"

85 Gallon

Description

070556 Mobile

070555 Mobile with Sump

070557 Stationary

- 16** • Dimensions: 46" x 24" x 25"

90 Gallon

Description

070559 Mobile

070558 Mobile with Sump

070560 Stationary

SPORTS TRAINER WHIRLPOOLS

(CON'T)

A FERNO-ILLE WHIRLPOOLS

- Designed specifically for the treatment and prevention of sport related injuries.
- Suggested use for treatment of arm, back, spine and leg injuries and are ideal for hamstring injuries.
- Includes a filling hose with faucet attachment.
- Dimensions: 48”L x 25”W x 28”D

Description	
071254	110 Gallon; Mobile
071255	110 Gallon; Stationary
071256	135 Gallon; Mobile
071257	135 Gallon; Stationary



A Features a second motor that powers the drain pump, completely emptying the tank.



B

SLANT BACK WHIRLPOOL

B SB SERIES

- Dimensions: 76” x 40” x 21”

100 Gallon	Description
070564	Mobile
070565	Stationary



C

Valves constantly adjust the outlet water to maintain the desired temperature.

WHIRLPOOL ACCESSORIES

C THERMOSTATIC MIXING VALVE

Whitehall	Description
070028	15 gpm
070029	25 gpm

Ferno-ille	Description
070069	10 gpm
070070	25 gpm

D



Adjustable flat seat is supported between “ladder” brackets hanging from tank rim.

D ADJUSTABLE SUSPENSION SEAT

Whitehall	Description
070032	For 20” Wide Tanks
070033	For 24” Wide Tanks
070034	For Low Tanks

Ferno-ille	Description
070073	For 22” Wide Tanks
070074	For 25” Wide Tanks

E



E COMBINATION TABLE WITH SEAT

Description
070086 For 20” Wide Tank; Whitehall
070087 For 24” Wide Tank; Whitehall

F

Seat Wraps around tank supported by brackets gripping the rim.



070079

Sturdy tubular chrome-plated steel frame facilitates treatment of the lower extremities.



070078

Sturdy tubular chrome-plated steel frame facilitates treatment of the lower extremities and features 3" ball-bearing casters, with two locking brakes.

G

070081



070080

**H**

Adding a small amount of this concentrated anti-foam emulsion before or during treatment controls excess foam. Non-toxic and non-irritating.

J

Slow-release chlorine that controls auto-and cross-contamination in treatment water. Approximately use 1 gram per gallon of water.

**I**

This concentrated, multi-purpose germicidal detergent and deodorant is effective against a broad spectrum of pathogenic bacteria, fungi and viruses.

**K**

Designed to be used with mobile and stationary extremity whirlpools, including High and Long and Low Whirlpools.

WHIRLPOOL ACCESSORIES (CON'T)

F TANK TOP SEAT

Whitehall	Description
070036	For 20" Wide Tank
070037	For 24" Wide Tank

Ferno-Illie	Description
070077	For 25" Wide Tank

G WHIRLPOOL CHAIRS

Low Whirlpool Chairs	Description
070078	Mobile; 33" to 43" Seat Height
070079	Stationary; 21" to 31.5" Seat Height

High Whirlpool Chairs	Description
070080	Mobile; 33" to 43" Seat Height
070081	Stationary; 28.75" to 38.75" Seat Height

WHIRLPOOL SUPPLIES

H D-FOAMING AGENT

- Compatible with all hydrotherapy additives.

Description
024017 One Gallon Bottle

I SANIZENE™

Description
024006 One Gallon Bottle
024007Z Case of (4) One Gallon Bottles

J CHLORAZENE® HYDROSOL SYSTEM

- Non-staining and non-foaming.

Description
024000 15 grams; (48) Packets/Box
024001 30 grams; (36) Packets/Box
024002 50 grams; (24) Packets/Box
024003 80 grams; (24) Packets/Box
024004 320 grams; (12) Bottles/Box
024005 5lb. Bottle

WHIRLPOOL PATIENT LIFTS

K MOBILE CHAIR LIFT



- Great as a positioning device to place your patient at a convenient working height for treatment.
- Rechargeable battery provides up to 40 continuous hours of use on a single charge.
- Hand-held controls allow the therapist to raise and lower the patient within the 32" range.
- Removable pads facilitate cleaning and disinfecting.
- Durable patient restraints and a handle provide security for the patient and move out of the way for patient access.
- Product weight: 225 lbs.

Length	Width	Height
070110 38.5"	34"	56"-88"

HYDRO THERAPY

WHIRLPOOL PATIENT LIFTS (CON'T)

A MOBILE STRETCHER LIFT

330lbs.
WEIGHT
CAPACITY

- Hand-held controls allow the therapist to raise and lower the patient within the 32" range.
- Rechargeable battery provides up to 40 continuous hours of use on a single charge.
- Removable pads facilitate cleaning and disinfecting.
- Durable patient restraints and a handle provide security for the patient and move out of the way for patient access.
- Product weight: 265 lbs.

	Length	Width	Height
070111	38.5"	78.75"	56"-88"

B AQUATIC ACCESS™ POOL LIFT

- Installs in minutes. Lifts 400lbs max.
- Portable, lifts out and rolls away with no tools
- Weighs 82lbs with seat, 57lbs without.
- Overall height in raised position 6'11"
- Compact, light weight, easy to remove and replace with no readjustments or tools required
- No underwater installation. Bottom of cylinder touches pool wall
- Protective vinyl-liner wall pad available.

	Description
070113	In Ground Pool Lift

AQUATIC THERAPY

C FERNO® AQUAGAITER™ FITNESS

- An underwater treadmill featuring a simple, chemical-resistant construction, utilizing a hydraulic direct-drive operating system.
- Concentrating on the natural buoyancy of the water to reduce joint impact, but with the resistance of the water against the body's movement, a safe aerobic workout is achieved.
- Suited for middle-aged and older adults, it is also great for athletes, the over weight, arthritic and the sports injured.
- Gentle on bones and heart.
- Hydrostatic pressure. Higher intensity at a lower heart rate.
- Running surface: 22" W x 58" L.
- Water depth: 2' to 5'.
- Product weight: 160 lbs.

	Length	Width	Height
070048	72"	34"	68"

A



Designed to be used with Hubbard tanks, burn tanks, and a wide variety of supine fixed and height-adjustable bathing and hydrotherapy units.

B



Hydrostatic pressure. Higher intensity at a lower heart rate.

E



HYDRO
THERAPY